

Lite Ryder

Lightweight Folding Powerchair

RC-Inn-N5909

.....



ROTHCARE



Features

.....

- Super Lightweight carbon fibre frame
- Folds in one easy step for transportation
- Comfortable padded seat & backrest
- Lightweight removable Lithium battery
- Compact easy to use joystick
- Complete chair weighs less than 13kg with the battery
- USB fitting in Joystick & in removable battery pack
- Battery pack is fitted with a charging port for offboard charging

Tech Specs

.....

Overall width	53cm
Overall length	90cm
Overall height	86cm
Seat width	43cm
Seat depth	40cm
Seat height from ground	46cm
Backrest height	40cm
Speed	4.5kph
Turning radius	90cm
Gradient	3 degrees
Range	10km
Battery	24V / 10Ah Lithium-ion

Motor	DC 24V, 150W brushless DC
Controller	Shanghai Micon
Brake	Electromagnetic
Front wheel	6"
Rear wheel	10.6"
Weight capacity	100kg
Weight w/out battery	11.85kg
Weight with battery	13kg

Local Provider: Emprise Mobility

Perth: 5 Sobek Pass | Mandurah: 14/8 Magenta Terrace
1300 761 196 | www.emprise.com.au

