

The Rea Azalea is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the seat position to be altered, ensuring maximum comfort at all times. This combined with good driving characteristics makes Azalea the ideal choice.







Rea Azalea

The Rea Azalea has seat tilt and a 30° backrest recline. Both tilt and recline can be easily adjusted.

Invacare's Dual Stability System (DSS) helps to prevent Decubitus thanks to its non-friction body positioning, which includes a biomechanical positioning backrest hinge that follows the user's natural back movement when adjusting the backrest angle.

Features



Flex 3, backrest plate. Apart from being height adjustable it is now width adjustable as well.



Option for electric adjustment of tilt and recline via a hand control that can be used by the attendant or user.

Tech Specs



Seat width 39 - 54cm



Total

Total width seat width + 23cm



Legrest angle Adjustable



Tilt Range 25°



Max. user weight



Transport weight 16.5kg



Seat height

45cm



77cm

lius Max. safe slope



7⁰

Cras App Fulf



Crash Tesst Approved Fulfills the requirements of ISO 7176-19



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com Thailand: Phone +66 0 2108 8555 or visit www.invacare.com

