

Mattress Selector for Adjustable Beds

Mattresses	Mattress style	Detail	Length	Sizes	SWL
I-Care IC15	Comfort mattress	Firm, memory foam, fabric cover	203cm	LS, KS, D, Q, SQ	200kg
I-Care IC20	Comfort mattress	Medium, memory foam, fabric cover	203cm	LS, KS, D, Q, SQ	200kg
I-Care IC25	Comfort mattress	Soft, memory foam, fabric cover	203cm	LS, KS, D, Q, SQ	200kg
I-Care M1	Medium Risk Bariatric Mattress	Castellated foam, firm edges, anti-shear cover	203cm	LS, KS, D, Q, SQ	300kg
I-Care M2	High Risk Medium Mattress	Castellated foam, firm edges, anti-shear cover	203cm	LS, KS, D, Q, SQ	215kg
I-Care M3	High Risk with Active Float Overlay	Castellated foam, overlay, firm edges, anti-shear cover	203cm	LS, KS, D, Q, SQ	215kg
Oska Fika	Very High Risk Soft Mattress	Multi-layer foam with heel relief, overlay, anti-shear cover	200cm	LS, KS	159kg
Oska Aria+	Very High Risk Foam / Air Combo	Foam in air cells, overlay, firm edges, anti-shear cover. Static or alternating.	200cm	LS, KS	250kg
Funke Hyper Foam 2	Very High Risk Mattress	Multi-layer foam with soft head and heel areas, firm edges, anti-shear cover	196cm	LS, KS	150kg
Funke Maxx 250	Very High Risk Mattress	Multi-layer foam with soft head and heel areas, firm edges, anti-shear cover	196cm	LS, KS	250kg
Funke Hyper Air	Very High Risk Foam / Air Combo	Foam in air cells, overlay, firm edges, anti-shear cover.	196cm	LS, KS	250kg
Funke ComfyCare	High Risk Soft Mattress	Multi-layer foam with soft head and heel areas, firm edges, fabric cover	203cm	LS, KS	150kg
Funke ComfyCare	High Risk Firm Mattress	Multi-layer foam with soft head and heel areas, firm edges, fabric cover	203cm	LS, KS	250kg
CairMax Duo	Very High Risk Foam / Air Combo	Foam in air cells, overlay, firm edges, anti-shear cover. Static or alternating.	200cm	LS, KS	230kg
Softform Premier	High Risk Medium Mattress	Castellated foam, firm edges, anti-shear cover	197cm	LS, KS	247kg
Softform MaxiGlide	High Risk Medium Mattress	Castellated foam, anti-shear layers, firm edges, anti-shear cover	197cm	LS, KS	247kg

Note: Adjustable beds are available in different lengths – check bed size to ensure the mattress is the right length.

All Icare beds are designed for mattresses 203cm long. Shorter mattresses will leave a gap at the head end of the bed but can still be selected.